

Erhöhung der Blutsicherheit: Neue UVC-basierte Technologie zur Pathogeninaktivierung von Thrombozytenkonzentraten zugelassen

1. Stramer SL, Hollinger FB, Katz LM, Kleinman S, Metzel PS, Gregory KR, Dodd RY: Emerging infectious disease agents and their potential threat to transfusion safety. *Transfusion* 2009;49 Suppl 2:1S-29S.
2. Orru S, Oberle D, Heiden M, Muller S, Krut O, Funk MB: Analysis of Transfusion-Transmitted Bacterial Infections according to German Hemovigilance Data (2011-2020). *Transfus Med Hemother* 2023;50:144-153.
3. Pitman JP, Payrat JM, Park MS, Liu K, Corash L, Benjamin RJ: Longitudinal analysis of annual national hemovigilance data to assess pathogen reduced platelet transfusion trends during conversion to routine universal clinical use and 7-day storage. *Transfusion* 2023;63:711-723.
4. Rasongles P, Angelini-Tibert MF, Simon P, Currie C, Isola H, Kientz D, Slaedts M, Jacquet M, Sundin D, Lin L, Corash L, Cazenave JP: Transfusion of platelet components prepared with photochemical pathogen inactivation treatment during a Chikungunya virus epidemic in Ile de La Reunion. *Transfusion* 2009;49:1083-1091.
5. Pierelli L, Vacca M, Zini G, Maresca M, Menichella G, Santinelli S, Stigliano MA, Mandarello G, Gasbarri R, Vaglio S: Emergency response of four transfusion centers during the last Chikungunya outbreak in Italy. *Transfusion* 2018;58:3027-3030.
6. Seltsam A, Mueller TH: UVC irradiation for pathogen reduction of platelet concentrates and plasma. *Transfus Med Hemother* 2011;38:43-54.
7. Faddy HM, Fryk JJ, Prow NA, Watterson D, Young PR, Hall RA, Tolksdorf F, Sumian C, Gravemann U, Seltsam A, Marks DC: Inactivation of dengue, chikungunya, and Ross River viruses in platelet concentrates after treatment with ultraviolet C light. *Transfusion* 2016;56:1548-1555.
8. Praditya D, Friesland M, Gravemann U, Handke W, Todt D, Behrendt P, Muller TH, Steinmann E, Seltsam A: Hepatitis E virus is effectively inactivated in platelet concentrates by ultraviolet C light. *Vox Sang* 2020;115:555-561.
9. Gravemann U, Handke W, Muller TH, Seltsam A: Bacterial inactivation of platelet concentrates with the THERAFLEX UV-Platelets pathogen inactivation system. *Transfusion* 2019;59:1324-1332.
10. Pohler P, Muller M, Winkler C, Schaudien D, Sewald K, Muller TH, Seltsam A: Pathogen reduction by ultraviolet C light effectively inactivates human white blood cells in platelet products. *Transfusion* 2015;55:337-347.
11. Brixner V, Bug G, Pohler P, Kramer D, Metzner B, Voss A, Casper J, Ritter U, Klein S, Alakel N, Peceny R, Derigs HG, Stegelmann F, Wolf M, Schrezenmeier H, Thiele T, Seifried E, Kapels HH, Doscher A, Petershofen EK, Muller TH, Seltsam A: Efficacy of UVC-treated, pathogen-reduced platelets versus untreated platelets: a randomized controlled non-inferiority trial. *Haematologica* 2021;106:1086-1096.
12. Mohr H, Gravemann U, Muller TH: Inactivation of pathogens in single units of therapeutic fresh plasma by irradiation with ultraviolet light. *Transfusion* 2009;49:2144-2151.
13. Handke W, Gravemann U, Muller TH, Wagner FF, Schulze TJ, Seltsam A: New ultraviolet C light-based method for pathogen inactivation of red blood cell units. *Transfusion* 2022;62:2314-2323.
14. Gregoire Y, Delage G, Custer B, Rochette S, Renaud C, Lewin A, Germain M: Cost-effectiveness of pathogen reduction technology for plasma and platelets in Quebec: A focus on potential emerging pathogens. *Transfusion* 2022;62:1208-1217.